

# FIRST RESPONSE<sup>®</sup>

## HCV Card Test

For the detection of antibodies to Hepatitis C virus in human serum / plasma / whole blood

### One-step Hepatitis C Card Test for Whole blood / Serum / Plasma

#### FEATURES

- ▶ Clear background.
- ▶ Self contained test procedure.
- ▶ Specifically engineered for whole blood specimen, serum & plasma.
- ▶ Storage temperature 4°C - 30°C.
- ▶ Detects HCV Subtypes like (1, 1a, 1b, 2, 2b, 3, 3a, 4, 5a and 6).
- ▶ Long Shelf life up to 24 months.

#### BENEFITS

- ▶ Unmistakable results at 15-20 minutes.
- ▶ Simple and easy to use, does not require advanced medical training.
- ▶ Cost effective, can be used effectively in resource-limited conditions.
- ▶ Accurate and reliable.
- ▶ Does not require cold chain & convenient for the small volume users.
- ▶ Ideal for laboratory / field conditions.



#### RESULTS INTERPRETATION

##### NEGATIVE




##### POSITIVE



##### INVALID RESULTS



#### TEST PROCEDURE



#### KIT CONTENTS

- Test device with desiccant
- Assay buffer bottle
- Specimen transfer device
- Sterile lancet
- Alcohol swab
- Instructions for use

#### EXTERNAL LABORATORY EVALUATION REPORTS

Place of Evaluation	Year of Testing	Sensitivity	Specificity
Zimbabwe	2019	100%	100%
NRL Melbourne, Australia	2020	100%	99.7%

#### ORDERING INFORMATION

info@premiermedcorp.com

Type : Bulk Packing

Catalogue No.	Pack Size
PI03FRC25	25 Tests / Kit
PI03FRC50	50 Tests / Kit
PI03FRC100	100 Tests / Kit
PI03FRC25CE	25 Tests / Kit
PI03FRC30CE	30 Tests / Kit



ISO 13485 & EN ISO 13485  
Certified Company



Manufactured by

**Premier Medical Corporation Private Limited**

Manufacturing Site and Registered Office : A1-302, GIDC, Sarigam 396155, Dist. Valsad, Gujarat, INDIA

• Customer support email : info@premiermedcorp.com • Tel.: +91 260 2780112/113,

• www.premiermedcorp.com