

FIRST RESPONSE[®]

Syphilis Anti-TP Card Test

For the detection of antibodies to *Treponema pallidum* in human whole blood / serum / plasma

One-step Anti-TP card test for Whole blood / Serum / Plasma



FEATURES

- ▶ Clear background
- ▶ Self contained test procedure
- ▶ Specifically engineered for whole blood, serum & plasma.
- ▶ 100% sensitivity and 100% specificity
- ▶ Storage temperature 4°C - 30°C
- ▶ Detects antibodies of all classes specific to *Treponema pallidum*
- ▶ Long Shelf life upto 24 months



BENEFITS


- ▶ Unmistakable results at 20 minutes.
- ▶ Simple and easy to use, does not require advanced medical training.
- ▶ Cost effective, can be used effectively in resource limited conditions.
- ▶ Accurate and reliable
- ▶ Does not require cold chain & convenient for small volume users.
- ▶ Ideal for laboratory / field conditions.




FIRST RESPONSE SYPHILIS



TEST PROCEDURE



RESULTS INTERPRETATION



KIT CONTENTS

- Test device with desiccant
- Sterile lancet
- Assay buffer bottle
- Alcohol swab
- Specimen transfer device
- Instructions for use

EXTERNAL LABORATORY EVALUATION REPORTS

Place of Evaluation	Year of Testing	Sensitivity	Specificity
Ministry of Health and Child Care, Zimbabwe	2019	100%	100%
Virus Research Institute, Uganda	2018	100%	100%
National Serology Reference Laboratory (NRL), Australia	2018	99.6%	100%

ORDERING INFORMATION

Info@premiermedcorp.com

Type : Bulk Packing

Catalogue No.	Pack Size
PI08FRC25	25 Tests / Kit
PI08FRC25-SA	25 Tests / Kit
PI08FRC30	30 Tests / Kit
PI08FRC50	50 Tests / Kit
PI08FRC100	100 Tests / Kit
PI08FRC25CE	25 Tests / Kit
PI08FRC30CE	30 Tests / Kit



ISO 13485 & EN ISO 13485
Certified Company

Part No. : PI08-BRO-001, Rev. AF, Date: 2025-04-01



Manufactured by

Premier Medical Corporation Private Limited

Mfg. Unit : A1-302, GIDC, Sarigam 396155. Dist. Valsad, Gujarat, INDIA

• Customer support email : info@premiermedcorp.com • Tel.: +91 260 2780112/113, • www.premiermedcorp.com