



**Malaria Ag. *P. falciparum* (HRP2/pLDH)
Combo Card Test (Whole Blood)**

The Clear Choice...

One Step Malaria Ag. *P. f.* (HRP2/pLDH) Combo


Detects Malaria *P. falciparum* specific to pLDH and HRP2 infection on single rapid card test (Whole blood)

FEATURES


- ▶ Detects *P. falciparum* specific to HRP2 & pLDH
- ▶ Self contained test procedure
- ▶ Detects as low as 200 parasites/μl.
- ▶ Test at any time of the day
- ▶ Long shelf life upto 24 Months.
- ▶ Storage temperature 1°C - 40°C.

BENEFITS

- ▶ Unmistakable results at 20 minutes.
- ▶ Differential diagnosis between *P. falciparum* specific to pLDH & HRP2
- ▶ Ideal for field testing
- ▶ Simple and easy to use.
- ▶ Cost effective
- ▶ Ideal in resource-limited condition
- ▶ Ideal for detecting low parasitemia levels.
- ▶ Allows immediate testing of specimen.



TEST PROCEDURE



KIT CONTENTS


- Test device with desiccant
- Sterile lancet
- Assay buffer bottle
- Alcohol swab
- Specimen transfer device
- Instructions for use

EXTERNAL EVALUATION DATA

Place of Evaluation	Year	Sensitivity 200 parasites/μl	Specificity
National Reference Laboratory	2024	100%	100%

RESULTS INTERPRETATION


Negative



Positive *P. f.* (pLDH)




Positive *P. f.* (HRP2)



Positive *P. f.* (pLDH & HRP2)



Invalid Results



ORDERING INFORMATION

info@premiermedcorp.com

Type : Bulk Packing

Catalogue No.	Pack Size
PI25FRC05	05 Tests / Kit
PI25FRC10	10 Tests / Kit
PI25FRC25	25 Tests / Kit
PI25FRC30	30 Tests / Kit

Type : Single Packing

PI25FRC10s	10 Single tests / Kit
PI25FRC25s	25 Single tests / Kit



**ISO 13485 & EN ISO 13485
Certified Company**

Part No. PI25-BRO-001, Rev. AA, Date:2025-07-10

Manufactured by



Premier Medical Corporation Private Limited

A1-302 and 3704-05, GIDC Sarigam, Valsad, Gujarat (India) - 396155

• Customer support email : info@premiermedcorp.com • Tel.: +91 260 2780112/113,
• www.premiermedcorp.com