

Malaria Ag. *P.falciparum* (HRP2/pLDH) Card Test (Whole Blood)

The Clear Choice...

One Step Malaria Ag. *P.f.* (HRP2/pLDH)

Detects Malaria *P. falciparum* specific to HRP2 and/or pLDH infection on single rapid card test (Whole blood)




FEATURES

- ▶ Detects *P. falciparum* specific to HRP2 and/or pLDH
- ▶ Self contained test procedure
- ▶ Detects as low as 200 parasites/μl.
- ▶ Test at any time of the day
- ▶ Long shelf life upto 24 Months.
- ▶ Storage temperature 1°C - 40°C.




BENEFITS

- ▶ Unmistakable results at 20 minutes.
- ▶ Ideal for field testing
- ▶ Simple and easy to use.
- ▶ Cost effective
- ▶ Ideal in resource-limited condition
- ▶ Ideal for detecting low parasitemia levels.
- ▶ Allows immediate testing of specimen.



TEST PROCEDURE



KIT CONTENTS

- Test device with desiccant
- Sterile lancet
- Assay buffer bottle
- Alcohol swab
- Specimen transfer device
- Instructions for use

EXTERNAL EVALUATION DATA


Place of Evaluation	Year	Sensitivity 200 parasites/μl	Specificity
National Reference Laboratory	2024	100%	100%

RESULTS INTERPRETATION

Negative



P.f. (HRP2 and/or pLDH) Positive



Invalid Results



ORDERING INFORMATION

info@premiermedcorp.com

Type : Bulk Packing

Catalogue No.	Pack Size
PI26FRC05	05 Tests / Kit
PI26FRC10	10 Tests / Kit
PI26FRC25	25 Tests / Kit
PI26FRC30	30 Tests / Kit

Type : Single Packing

PI26FRC10s	10 Single tests / Kit
PI26FRC25s	25 Single tests / Kit



ISO 13485 & EN ISO 13485
Certified Company

Part No. PI26-BRO-001, Rev. AA, Date:2025-09-08



Manufactured by

Premier Medical Corporation Private Limited

A1-302 and 3704-05, GIDC Sarigam, Valsad, Gujarat (India) - 396155

• Customer support email : info@premiermedcorp.com • Tel.: +91 260 2780112/113,
• www.premiermedcorp.com