



# FIRST RESPONSE

Malaria Ag. *P.f.* (HRP2) / *P.f.* (pLDH) / *P.v.* (pLDH) Test (Whole Blood)

## The Clear Choice...

### Malaria Ag. *P.f.* (HRP2) | *P.f.* (pLDH) | *P.v.* (pLDH) Combo


Detects Malaria Ag. *P.falciparum* specific to pLDH and HRP2 and *P.vivax* specific to pLDH in human whole blood

#### FEATURES


- ▶ Detects Malaria Ag. *P.falciparum* specific to pLDH and HRP2 and *P.vivax* specific to pLDH in human whole blood
- ▶ Self contained test procedure
- ▶ Detects as low as 200 parasites/μl.
- ▶ Test at any time of the day
- ▶ Long shelf life upto 24 Months.
- ▶ Storage temperature 1°C - 40°C.

#### BENEFITS

- ▶ Unmistakable results at 20 minutes. Differential diagnosis between *P.falciparum* specific to pLDH & HRP2 & *P.vivax* specific to pLDH
- ▶ Ideal for field testing
- ▶ Simple and easy to use.
- ▶ Cost effective
- ▶ Ideal in resource-limited condition
- ▶ Ideal for detecting low parasitemia levels.
- ▶ Allows immediate testing of specimen.



#### TEST PROCEDURE



#### KIT CONTENTS

- Test device with desiccant
- Sterile lancet
- Assay buffer bottle
- Alcohol swab
- Specimen transfer device
- Instructions for use

#### ORDERING INFORMATION

info@premiermedcorp.com

Type : Bulk Packing

Catalogue No.	Pack Size
PI29FRC05	05 Tests / Kit
PI29FRC10	10 Tests / Kit
PI29FRC25	25 Tests / Kit
PI29FRC30	30 Tests / Kit

Type : Single Packing

PI29FRC10s	10 Single tests / Kit
PI29FRC25s	25 Single tests / Kit

#### EXTERNAL EVALUATION DATA

Place of Evaluation	Year	Sensitivity 200 parasites/μl	Specificity
National Reference Laboratory	2025	100%	100%

#### RESULTS INTERPRETATION

Negative	<i>P.f.</i> (HRP2 & pLDH) Positive	<i>P.f.</i> (HRP2) Positive	<i>P.f.</i> (pLDH) Positive	<i>P.v.</i> Positive	<i>P.f.</i> (HRP2), <i>P.f.</i> (pLDH) & <i>P.v.</i> (pLDH) Positive	<i>P.f.</i> (pLDH) & <i>P.v.</i> (pLDH) Positive	<i>P.f.</i> (HRP2) & <i>P.v.</i> (pLDH) Positive	Invalid Results

Manufactured by

**Premier Medical Corporation Private Limited**

A1-302 and 3704-05, GIDC Sarigam, Valsad, Gujarat (India) - 396155

• Customer support email : info@premiermedcorp.com • Tel.: +91 260 2780112/113,  
• www.premiermedcorp.com

An ISO 13485 & EN ISO 13485 Certified Company

Part No. PI29-BRO-001, Rev. AA, Date:2026-06-05