



Premier Medical Corporation's

Triple Screen Program

HIV + Syphilis + Hepatitis B Screening Solution



Features & Benefits

One Program. Three Critical Tests.

An integrated Triple Screen Program delivering laboratory-comparable

HIV, Syphilis, and HBsAg results in a single antenatal visit—enabling early diagnosis, timely treatment, and birth-dose vaccination to prevent mother-to-child transmission.”

HIV + Syphilis Combo

- 100% sensitivity & specificity
- Detects gp41 & gp36 of HIV
- Detects P47, P45, P17 & P15 of Syphilis
- Clear differentiation between HIV & Syphilis infection
- Result in 15 minutes

HBsAg Test

- 100% sensitivity & specificity
- Detects all subtypes of HBsAg
- Result in 20 minutes
- Detect HBV early
- Prevent transmission at birth.

Common Features

- Specimen Type: Whole Blood/Serum/Plasma
- Easy to use, no advanced training required
- High accuracy & reliable results
- No cold chain required
- Suitable for lab & field conditions
- One Visit. Three Tests. Zero Transmission.

Why Triple Screening?

- Early detection of HIV, Syphilis and HBV
- Saves time with combined testing
- Prevents mother-to-child transmission
- Cost-effective for large-scale antenatal programs
- Ideal for NGOs, camps & hospitals
- Supports national elimination goals

Digital Support Tools



Assaya Digital Reader

- Enables digital reading of RDT results
- High sensitivity with GPS tracking
- Supports real-time infection mapping



Scan to Visit
Assaya Digital Reader



Learning Management Software (LMS)

- Online training for healthcare users
- Helps in correct test usage & interpretation
- Direct support from Premier Medical team



Scan to Visit
Learning Management Software



Premier Medical Corporation Private Limited

A1-302 and 3704-05, GIDC, Sarigam 396155.
Dist. Valsad, Gujarat, INDIA

• Customer support email : info@premiermedcorp.com
• Tel.: +91 260 2780112/113, • www.premiermedcorp.com