



**1** **FIRST**<sup>®</sup>  
**RESPONSE**

## HBsAg Card Test

For the detection of Hepatitis B surface antigen in human whole blood / serum / plasma

### One-step Hepatitis B Card Test for Whole blood / Serum / Plasma

#### FEATURES

- ▶ Clear background.
- ▶ Self contained test procedure.
- ▶ Specifically engineered for whole blood specimen, serum & plasma.
- ▶ 100% sensitivity and 100% specificity.
- ▶ Storage temperature 4°C - 30°C.
- ▶ Detects all subtypes of HBsAg.
- ▶ Long shelf life up to 24 months.

#### BENEFITS

- ▶ Unmistakable results at 20 minutes.
- ▶ Simple and easy to use, does not require advanced medical training.
- ▶ Cost effective, can be used effectively in resource limited conditions.
- ▶ Accurate and reliable.
- ▶ Does not require cold chain & convenient for small volume users.
- ▶ Ideal for laboratory / field conditions.




#### RESULTS INTERPRETATION

##### NEGATIVE




##### POSITIVE



##### INVALID RESULTS



#### TEST PROCEDURE



#### KIT CONTENTS

- Test device with desiccant
- Assay buffer bottle
- Specimen transfer device
- Sterile lancet
- Alcohol swab
- Instructions for use

#### EXTERNAL LABORATORY EVALUATION REPORTS

Place of Evaluation	Year of Testing	Sensitivity	Specificity
Biosynex, France	2018	100%	100%
Virus Research Institute, Uganda	2018	99.1%	100%

#### ORDERING INFORMATION

info@premiermedcorp.com

Type : Bulk Packing

Catalogue No.	Pack Size
PI10FRC25	25 Tests / Kit
PI10FRC30	30 Tests / Kit
PI10FRC25CE	25 Tests / Kit
PI10FRC30CE	30 Tests / Kit



ISO 13485 & EN ISO 13485  
Certified Company

Part No. : PI10-BRO-001, Rev. AA, Date: 2021-01-09



Manufactured by

**Premier Medical Corporation Private Limited**

A1-302, GIDC, Sarigam 396155, Dist. Valsad, Gujarat, INDIA

• Customer support email : info@premiermedcorp.com • Tel.: +91 260 2780112/113,  
• www.premiermedcorp.com