


HCV Card Test

For the detection of antibodies to Hepatitis C virus in human serum / plasma / whole blood



One-step Hepatitis C Card Test for Whole blood / Serum / Plasma

FEATURES


- ▶ Clear background.
- ▶ Self contained test procedure.
- ▶ Specifically engineered for whole blood specimen, serum & plasma.
- ▶ 100% sensitivity and 100% specificity.
- ▶ Storage temperature 4°C - 30°C.
- ▶ Detects all subtypes of HCV.
- ▶ Long Shelf life up to 24 months.

BENEFITS

- ▶ Unmistakable results at 15-20 minutes.
- ▶ Simple and easy to use, does not require advanced medical training.
- ▶ Cost effective, can be used effectively in resource-limited conditions.
- ▶ Accurate and reliable.
- ▶ Does not require cold chain & convenient for the small volume users.
- ▶ Ideal for laboratory / field conditions.



TEST PROCEDURE



KIT CONTENTS

- Test device with desiccant
- Assay buffer bottle
- Specimen transfer device
- Sterile lancet
- Alcohol swab
- Instructions for use

EXTERNAL LABORATORY EVALUATION REPORTS

Place of Evaluation	Year of Testing	Sensitivity	Specificity
Zimbabwe	2019	100%	100%
NRL Melbourne, Australia	2020	100%	99.7%

RESULTS INTERPRETATION


NEGATIVE



POSITIVE



INVALID RESULTS



ORDERING INFORMATION

info@premiermedcorp.com
Type : Bulk Packing

Catalogue No.	Pack Size
PI03FRC25	25 Tests / Kit
PI03FRC30	30 Tests / Kit
PI03FRC25CE	25 Tests / Kit
PI03FRC30CE	30 Tests / Kit