



Prequalified



# First Response Malaria Antigen *P. falciparum* (HRP2) Card Test


## The Clear Choice...

### One Step MALARIA Antigen *P. falciparum* (HRP2) Detection Card Test

Detects *P. falciparum* infection on rapid card test (Whole blood)

#### FEATURES

- ▶ Clear background
- ▶ Self contained test procedure
- ▶ 100% sensitivity and 100% specificity
- ▶ Storage temperature 1°C - 40°C
- ▶ Detection of *P. falciparum* in human blood, accurate at parasite concentrations of 200 parasites/μl.
- ▶ High tech membrane strip locks in visible test results at ambient conditions even months later.
- ▶ Long shelf life up to 24 Months




#### BENEFITS

- ▶ Unmistakable results at 20 minutes
- ▶ Simple and easy to use, ideal in resource-limited conditions or blood bank screening
- ▶ Accurate and reliable
- ▶ Ships as normal cargo, no requirement for cold chain, direct and indirect cost benefits
- ▶ Early detection of most virulent *Plasmodium falciparum*.
- ▶ Convenient for small volume users.


#### KIT CONTENTS

- Test device with desiccant
- Assay buffer bottle
- Specimen transfer device
- Sterile lancet
- Alcohol swab
- Instructions for use

#### TEST PROCEDURE



#### RESULTS INTERPRETATION



#### EXTERNAL EVALUATION DATA

Name of the Institute	Year of Testing	Specificity	Sensitivity	
			200 P/μl	2000 P/μl
WHO Evaluation Round 6	2014-2015	100%	91%	100%
National Public Health & Reference Lab, Ghana	2015	100%	100%	

#### ORDERING INFORMATION

Info@premiermedcorp.com	
Type : Bulk Packing	
Catalogue No.	Pack Size
PI13FRC05	05 Tests / Kit
PI13FRC25	25 Tests / Kit
PI13FRC05CE	05 Tests / Kit
PI13FRC25CE	25 Tests / Kit
Type : Single Packing	
PI13FRC10s	10 Single tests / Kit
PI13FRC25s	25 Single tests / Kit

